Table 1: 测量普朗克常数 h 及红限频率数据表

| 波长λ/nm | 365 | 405 | 436 | 546 | 577 | | |
|--------------|--------|--------|------------------|--------|--------|--|--|
| $U_s/{ m V}$ | -1.458 | -1.300 | -1.003 -1.035 | -0.434 | -0.313 | | |
| | | | | -0.456 | -0.336 | | |
| | -1.495 | -1.326 | -10.43 | -0.468 | -0.346 | | |
| | -1.502 | -1.335 | -1.052 | -0.477 | -0.365 | | |
| 平均值/V | -1.482 | -1.321 | -1.033 | -0.459 | -0.338 | | |




Table 2: 测量光电管伏安特性数据表

| 光阑/mm | 2 | | 4 | | 8 | |
|-------|------|---------------|------|---------------|------|---------------|
| 测量量 | U/V | $I/10^{-13}A$ | U/V | $I/10^{-12}A$ | U/V | $I/10^{-12}A$ |
| | -0.4 | 0.00 | -0.3 | 0.00 | -0.3 | 0.00 |
| | 1.1 | 577 | 1.2 | 180.2 | 1.2 | 383 |
| | 2.6 | 908 | 2.7 | 336 | 2.7 | 779 |
| | 4.1 | 1114 | 4.2 | 408 | 4.2 | 1129 |
| | 5.6 | 1293 | 5.7 | 479 | 5.7 | 1406 |
| | 7.1 | 1460 | 7.2 | 543 | 7.2 | 1623 |
| 数据 | 8.6 | 1620 | 8.7 | 602 | 8.7 | 1795 |
| | 10.1 | 1768 | 10.2 | 657 | 10.2 | 1948 |
| | 11.6 | 1889 | 11.7 | 707 | 11.7 | 2090 |
| | 13.1 | 1997 | 13.2 | 758 | 13.2 | 2270 |
| | 14.6 | 2150 | 14.7 | 799 | 14.7 | 2390 |
| | 16.1 | 2230 | 16.2 | 835 | 16.2 | 2520 |
| | 17.6 | 2310 | 17.7 | 863 | 17.7 | 2640 |
| | 19.1 | 2370 | 19.2 | 906 | 19.2 | 2720 |
| | 20.6 | 2470 | 20.7 | 932 | 20.7 | 2810 |
| | 22.1 | 2530 | 22.2 | 958 | 22.2 | 2870 |
| | 23.6 | 2580 | 23.7 | 981 | 23.7 | 2940 |
| | 25.1 | 2660 | 25.2 | 998 | 25.2 | 3030 |
| | 26.6 | 2700 | 26.7 | 1021 | 26.7 | 3080 |
| | 28.1 | 2750 | 28.2 | 1050 | 28.2 | 3170 |
| | 29.6 | 2800 | 29.7 | 1062 | 29.7 | 3230 |

Table 3: 验证光饱和电流与入射光强正比关系数据表

| 波长/nm | | 数据 | | |
|-------|---|--|-----------|-----------|
| 546 | 光阑/mm <i>I/</i> 10 ⁻¹² <i>A</i> | $\begin{array}{ c c c }\hline 2\\174.3\end{array}$ | $4\\669$ | 8 2660 |
| 436 | 光阑/mm <i>I/</i> 10 ⁻¹² <i>A</i> | 2 514 | 4 1879 | 8 6860 |